

TargetMol®
YOUR TARGET MOLECULES

Chemical Properties

C1CN(CCN1Cc2c[nH]c3ccccc23)CCN4CCCCC4

Description	ZEL-H16 is an agonist of the histamine H3 receptor.
-------------	---

	1mg	5mg	10mg
1 mM	2.7206 mL	13.6029 mL	27.2057 mL
5 mM	0.5441 mL	2.7206 mL	5.4411 mL
10 mM	0.2721 mL	1.3603 mL	2.7206 mL
50 mM	0.0544 mL	0.2721 mL	0.5441 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481