

VU0155056

Chemical Properties


CAS No. : 1130067-18-3

Formula: C₂₅H₂₆N₄O₂

Molecular Weight: 414.5

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	VU0155056 is a cell-permeable phospholipase D1 and D2 inhibitor.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4125 mL	12.0627 mL	24.1255 mL
5 mM	0.4825 mL	2.4125 mL	4.8251 mL
10 mM	0.2413 mL	1.2063 mL	2.4125 mL
50 mM	0.0483 mL	0.2413 mL	0.4825 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481