

U 19052

Chemical Properties


CAS No. : 112682-32-3

Formula: C₃₁H₄₂O₆S

Molecular Weight: 542.73

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	U 19052 is a Leukotriene D4 antagonist.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8425 mL	9.2127 mL	18.4254 mL
5 mM	0.3685 mL	1.8425 mL	3.6851 mL
10 mM	0.1843 mL	0.9213 mL	1.8425 mL
50 mM	0.0369 mL	0.1843 mL	0.3685 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481