


MDL-27088

Chemical Properties

CAS No. : 111581-85-2
Formula: C₂₅H₂₈N₂O₅
Molecular Weight: 436.5
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	MDL-27088 is an angiotensin-converting enzyme inhibitor.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.291 mL	11.4548 mL	22.9095 mL
5 mM	0.4582 mL	2.291 mL	4.5819 mL
10 mM	0.2291 mL	1.1455 mL	2.291 mL
50 mM	0.0458 mL	0.2291 mL	0.4582 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481