


Amphidinolide D

Chemical Properties

CAS No. : 110786-78-2
Formula: C₃₂H₅₀O₈
Molecular Weight: 562.73
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Amphidinolide D is a powerful activator of actomyosin ATPase enhances skeletal muscle contraction.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7771 mL	8.8853 mL	17.7705 mL
5 mM	0.3554 mL	1.7771 mL	3.5541 mL
10 mM	0.1777 mL	0.8885 mL	1.7771 mL
50 mM	0.0355 mL	0.1777 mL	0.3554 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481