


Pincaïnide

Chemical Properties

CAS No. : 83471-41-4
Formula: C₁₆H₂₄N₂O
Molecular Weight: 260.37
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Pincaïnide is a novel beta-aminoaniline that inhibits, in a dose-dependent manner, diseases caused by abnormal norepinephrine vascular smooth muscle function.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8407 mL	19.2034 mL	38.4069 mL
5 mM	0.7681 mL	3.8407 mL	7.6814 mL
10 mM	0.3841 mL	1.9203 mL	3.8407 mL
50 mM	0.0768 mL	0.3841 mL	0.7681 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Bourghardt Peebo B, et al. Expression of the focal adhesion protein PINCH in normal and alkali-injured corneas and the role of PMNs. Acta Ophthalmol Scand. 2007 Jun;85(4):395-400.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481