


Data Sheet (Cat.No.T71964)

Ravidasvir HCl

Chemical Properties

CAS No. : 1303533-81-4
Formula: C42H52Cl2N8O6
Molecular Weight: 835.828
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Ravidasvir, also known as PPI-668 and ASC16, is a second-generation, orally active, potent and selective HCV NS5A protein inhibitor.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1964 mL	5.9821 mL	11.9642 mL
5 mM	0.2393 mL	1.1964 mL	2.3928 mL
10 mM	0.1196 mL	0.5982 mL	1.1964 mL
50 mM	0.0239 mL	0.1196 mL	0.2393 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481