


PKC-9

Chemical Properties

CAS No. : 1132609-87-0
Formula: C₂₅H₂₅N₇
Molecular Weight: 423.51
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	PKC-9 is a PKC-zeta inhibitor 9.
-------------	----------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3612 mL	11.8061 mL	23.6122 mL
5 mM	0.4722 mL	2.3612 mL	4.7224 mL
10 mM	0.2361 mL	1.1806 mL	2.3612 mL
50 mM	0.0472 mL	0.2361 mL	0.4722 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481