

FXR agonist 4

Chemical Properties


CAS No. :

Formula: C₂₁H₂₈ClN₃O

Molecular Weight: 373.92

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	FXR agonist 4 with an EC ₅₀ value of 1.05 μM. FXR agonist 4 effectively improves hyperlipidemia, hepatic steatosis, insulin resistance and hepatic inflammation in DIO mice. FXR agonist 4 can be used for the research of non-alcoholic fatty liver disease (NAFLD) .
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6744 mL	13.3718 mL	26.7437 mL
5 mM	0.5349 mL	2.6744 mL	5.3487 mL
10 mM	0.2674 mL	1.3372 mL	2.6744 mL
50 mM	0.0535 mL	0.2674 mL	0.5349 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481