


METTL3-IN-1

Chemical Properties

CAS No. : 2648257-77-4
Formula: C₂₅H₂₄FN₅O₂
Molecular Weight: 445.49
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	METTL3-IN-1 is a polyheterocyclic compound, acts as METTL3 inhibitor .
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2447 mL	11.2236 mL	22.4472 mL
5 mM	0.4489 mL	2.2447 mL	4.4894 mL
10 mM	0.2245 mL	1.1224 mL	2.2447 mL
50 mM	0.0449 mL	0.2245 mL	0.4489 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481