


BAY-3153

Chemical Properties

CAS No. : 2771319-69-6
Formula: C₂₅H₂₉Cl₂N₃O₄
Molecular Weight: 506.42
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	BAY-3153 is a selective CCR1 (C-C motif chemokine receptor 1) antagonist (human IC ₅₀ =3 nM ; rat IC ₅₀ =11 nM ; mice IC ₅₀ =81 nM) .
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9746 mL	9.8732 mL	19.7465 mL
5 mM	0.3949 mL	1.9746 mL	3.9493 mL
10 mM	0.1975 mL	0.9873 mL	1.9746 mL
50 mM	0.0395 mL	0.1975 mL	0.3949 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481