

Peginterferon beta-1a

Chemical Properties

CAS No. : 1211327-92-2

Formula:

Molecular Weight:

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	Peginterferon beta-1a (Peginterferon β -1a), the inaugural pegylated interferon beta-1a derivative, promotes apoptosis in cancer cells and exhibits anti-tumor properties in nude mice models. It is utilized in the study of cancer and relapsing multiple sclerosis (RMS) [1].
-------------	--

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481