

Peginterferon alfa-2b

Chemical Properties

CAS No. : 215647-85-1
Formula:
Molecular Weight:
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	Peginterferon alfa-2b (Sch 54031), a recombinant alfa-2b interferon covalently bonded to polyethylene glycol (PEG), acts as an immunomodulator. This compound is effective in reducing viral DNA in HIV and has potential applications in melanoma research [1] [2].
-------------	--

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481