


TGFβR-IN-1

Chemical Properties

CAS No. : 2857845-36-2
Formula: C₈₄H₁₂₁ClFN₁₇O₂₃
Molecular Weight: 1791.41
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	TGFβR-IN-1 is a long-acting tumor-activated proagent of a TGFβR inhibitor.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.5582 mL	2.7911 mL	5.5822 mL
5 mM	0.1116 mL	0.5582 mL	1.1164 mL
10 mM	0.0558 mL	0.2791 mL	0.5582 mL
50 mM	0.0112 mL	0.0558 mL	0.1116 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481