


[18F]-Labeled L-dopa precursor

Chemical Properties

CAS No. : 1614253-58-5
Formula: C₃₅H₃₆N₂O₇
Molecular Weight: 596.67
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	[18F]-Labeled L-dopa precursor is a precursor for synthesis of 18F-labeled L-dopa[1].
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.676 mL	8.3798 mL	16.7597 mL
5 mM	0.3352 mL	1.676 mL	3.3519 mL
10 mM	0.1676 mL	0.838 mL	1.676 mL
50 mM	0.0335 mL	0.1676 mL	0.3352 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481