

Milpocitide

Chemical Properties

CAS No. : 2643306-81-2
Formula: C255H413N63O86S6
Molecular Weight: 5929.75
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	Milpocitide, a (293-333)-peptide fragment of the human low-density lipoprotein receptor (LDLR), specifically corresponds to the EGF-like domain 1 [1].
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.1686 mL	0.8432 mL	1.6864 mL
5 mM	0.0337 mL	0.1686 mL	0.3373 mL
10 mM	0.0169 mL	0.0843 mL	0.1686 mL
50 mM	0.0034 mL	0.0169 mL	0.0337 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481