Data Sheet (Cat.No.T74637)



Sonrotoclax

Chemical Properties


CAS No.: 2383086-06-2

Formula: C49H59N7O7S

Molecular Weight: 890.1

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description Sonrotoclax is a potent, orally active inhibitor of Bcl2, demonstrating effective cell-

killing activity against a range of lymphoma and leukemia cell lines [1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1235 mL	5.6173 mL	11.2347 mL
5 mM	0.2247 mL	1.1235 mL	2.2469 mL
10 mM	0.1123 mL	0.5617 m ⊡	1.1235 mL
50 mM	0.0225 mL	0.1123 mL	0.2247 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481

Page 1 of 1 www.targetmol.com