

Data Sheet (Cat.No.T74868)

Uracil-m7GpppAmpG ammonium

Chemical Properties

CAS No. :

Formula: C₃₁H₄₂N₁₂O₂₅P₄.xNH₃

Molecular Weight:

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	Uracil-m7GpppAmpG ammonium is a cap analog that can be used for mRNA synthesis.
-------------	---

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481