

Data Sheet (Cat.No.T74872)

ODN 105870

Chemical Properties

CAS No. : 1965342-40-8

Formula:

Molecular Weight:

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	ODN 105870, a G-modified inhibitory oligonucleotide (INH-ODN), selectively inhibits TLR7 and is utilized in the research of inflammatory immune responses [1].
-------------	--

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481