

WU-FA-01

Chemical Properties

CAS No. : 882429-53-0
Formula: C₃₄H₅₂O₉
Molecular Weight: 604.77
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	WU-FA-01, a hydrogenated derivative of WU-FA-00, functions as an antibacterial agent with robust activity against Gram-positive strains and possesses anti-inflammatory properties [1].
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6535 mL	8.2676 mL	16.5352 mL
5 mM	0.3307 mL	1.6535 mL	3.307 mL
10 mM	0.1654 mL	0.8268 mL	1.6535 mL
50 mM	0.0331 mL	0.1654 mL	0.3307 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481