

ATU027

Chemical Properties

CAS No. : 1415935-22-6

Formula:

Molecular Weight:

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	ATU027, a small interfering RNA (siRNA), silences protein kinase N3 (PKN3) expression in the vascular endothelium. It has been demonstrated to thwart local tumor invasion, as well as lymph node and pulmonary metastasis, in mouse cancer models.
-------------	---

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481