

EGF Receptor Substrate 2 (Phospho-Tyr5)

Chemical Properties

CAS No. :

Formula: C54H82N13O24

Molecular Weight: 1328.28

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	EGF Receptor Substrate 2 (Phospho-Tyr5) is a biologically active peptide derived from an autophosphorylation site (Tyr 992) of epidermal growth factor receptor (EGFR).
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7529 mL	3.7643 mL	7.5285 mL
5 mM	0.1506 mL	0.7529 mL	1.5057 mL
10 mM	0.0753 mL	0.3764 mL	0.7529 mL
50 mM	0.0151 mL	0.0753 mL	0.1506 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481