

(Asp28)-Exenatide

Chemical Properties

CAS No. :	1678417-24-7
Formula:	C184H281N49O61S
Molecular Weight:	4187.56
Appearance:	no data available
Storage:	keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year

Biological Description

Description	(Asp28)-Exenatide is a degradation product of exenatide . (Asp28)-Exenatide can be used as a GLP-1R agonist [1] .
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.2388 mL	1.194 mL	2.388 mL
5 mM	0.0478 mL	0.2388 mL	0.4776 mL
10 mM	0.0239 mL	0.1194 mL	0.2388 mL
50 mM	0.0048 mL	0.0239 mL	0.0478 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481