

Lebrikizumab

Chemical Properties

CAS No. : 953400-68-5
Formula:
Molecular Weight:
Appearance: no data available
Storage: Powder -20°C for 3 years
In solvent -80°C for 2 years

Biological Description

Description	Lebrikizumab, an IgG4 humanized monoclonal antibody, specifically targets and inhibits interleukin-13 (IL-13), thus impacting its function. This compound is utilized in asthma research [1].
-------------	---

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481