

Lenzilumab

Chemical Properties

CAS No. : 1229575-09-0
Formula:
Molecular Weight:
Appearance: no data available
Storage: Powder -20°C for 3 years
In solvent -80°C for 2 years

Biological Description

Description	Lenzilumab (KB 003) is a humanized monoclonal antibody designed to target CSF2/GM-CSF, with applications in COVID-19 and studies on chronic myelomonocytic leukemia (CMML) and juvenile myelomonocytic leukemia (JMML) [1] [2].
-------------	---

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481