

Befovacimab

Chemical Properties

CAS No. : 2156634-62-5

Formula:

Molecular Weight:

Appearance: no data available

Storage:

Powder	-20°C for 3 years
In solvent	-80°C for 2 years

Biological Description

Description	Befovacimab (BAY 1093884), a fully human monoclonal IgG2 antibody, binds to the tissue factor pathway inhibitor (TFPI). It is utilized in research for haemophilia A/B [1].
-------------	---

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481