

Iscalimab

Chemical Properties

CAS No. : 2031153-61-2

Formula:

Molecular Weight:

Appearance: no data available

Storage: Powder -20°C for 3 years
In solvent -80°C for 2 years

Biological Description

Description	Iscalimab (CFZ-533), a non-depleting IgG1 monoclonal antibody that targets CD40 with a dissociation constant (K _D) of 0.3 nM, is utilized in the research of Graves' hyperthyroidism and autoimmune diseases [1] [2] [3].
-------------	---

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481