

Regavirumab

Chemical Properties

CAS No. : 153101-26-9

Formula:

Molecular Weight:

Appearance: no data available

Storage: Powder -20°C for 3 years
In solvent -80°C for 2 years

Biological Description

Description	Regavirumab (MCA C23), a human monoclonal antibody, targets cytomegalovirus (CMV), human serum albumin (HSA), and amino acetic acid [1].
-------------	--

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481