

Data Sheet (Cat.No.T77273)

DSPE-PEG-Azide (MW 2000)

Chemical Properties

CAS No. :

Formula: 2000(Average)

Molecular Weight:

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

Description	DSPE-PEG-Azide (MW 2000) is an azide containing lipid that can be used to form micelles as nanoparticles for drug delivery [1] .
-------------	--

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481