

DSPE-alkyne

Chemical Properties


CAS No. :

Formula: C45H84NO10P

Molecular Weight: 830.12

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	DSPE-alkyne is a phospholipid [1] .
-------------	-------------------------------------

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2046 mL	6.0232 mL	12.0465 mL
5 mM	0.2409 mL	1.2046 mL	2.4093 mL
10 mM	0.1205 mL	0.6023 mL	1.2046 mL
50 mM	0.0241 mL	0.1205 mL	0.2409 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481