


TV 3279

Chemical Properties

CAS No. : 209394-29-6
Formula: C₁₆H₂₀N₂O₂
Molecular Weight: 272.34
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	TV 3279 is a novel ChE-MAI inhibitor , and the neuroprotective properties depend on their ability to induce the anti-apoptotic proteins PKC, Bcl-2, Bcl-x, and SOD, and to block the nuclear translocation of the pro-apoptotic enzyme glyceraldehyde phosphate dehydrogenase in PC-12 and neuroblastoma cells.
Targets(IC50)	BCL,PKC
In vivo	Chronic TV3279 (26 mg kg ⁻¹ for 21 days; rats) similarly inhibited 50% of striatal ChE activity but did not affect MAO activity or amine levels. Accordingly, chronic TV3279 abolished the ability of TCP to initiate L-dopa-induced hyperactivity. In contrast, chronic as well as acute treatments with TV3279 reduced spontaneous motor activity.[2]

Solubility Information

Solubility	DMSO: 2.73 mg/mL (10.01 mM) (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6719 mL	18.3594 mL	36.7188 mL
5 mM	0.7344 mL	3.6719 mL	7.3438 mL
10 mM	0.3672 mL	1.8359 mL	3.6719 mL
50 mM	0.0734 mL	0.3672 mL	0.7344 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Sterling J, et al. Novel dual inhibitors of AChE and MAO derived from hydroxy aminoindan and phenethylamine as potential treatment for Alzheimer's disease. J Med Chem. 2002;45(24):5260-5279.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481