


4-Methyl-6-Methoxycoumarin

Chemical Properties

CAS No. : 6295-35-8
Formula: C₁₁H₁₀O₃
Molecular Weight: 190.2
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	4-Methyl-6-Methoxycoumarin is a natural product isolated from <i>Eupatorium pauciflorum</i> .
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 25 mg/mL (131.44 mM), (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.2576 mL	26.2881 mL	52.5762 mL
5 mM	1.0515 mL	5.2576 mL	10.5152 mL
10 mM	0.5258 mL	2.6288 mL	5.2576 mL
50 mM	0.1052 mL	0.5258 mL	1.0515 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481