


(+)-Dehydroabietic acid

Chemical Properties

CAS No. :	1231-75-0
Formula:	C ₂₀ H ₂₈ O ₂
Molecular Weight:	300.4
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	(+)-Dehydroabietic acid is a natural compound isolated from Rosin acid.
Targets(IC ₅₀)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3289 mL	16.6445 mL	33.2889 mL
5 mM	0.6658 mL	3.3289 mL	6.6578 mL
10 mM	0.3329 mL	1.6644 mL	3.3289 mL
50 mM	0.0666 mL	0.3329 mL	0.6658 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481