

BCL6-IN-7

Chemical Properties


CAS No. : 2097518-46-0

Formula: C₁₈H₁₅ClN₆O

Molecular Weight: 366.8

Appearance: no data available

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	BCL6-IN-7 is a potent BCL6 protein-protein interaction(PPI) inhibitor[1].
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7263 mL	13.6314 mL	27.2628 mL
5 mM	0.5453 mL	2.7263 mL	5.4526 mL
10 mM	0.2726 mL	1.3631 mL	2.7263 mL
50 mM	0.0545 mL	0.2726 mL	0.5453 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481