


1,2-Distearoyl-sn-glycero-3-phosphorylcholine

Chemical Properties

CAS No. :	816-94-4
Formula:	C44H88NO8P
Molecular Weight:	790.15
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	1,2-Distearoyl-sn-glycero-3-phosphorylcholine (L-ALPHA-PHOSPHATIDYLCHOLINE, DISTEAROYL) is a cylindrical-shaped lipid used to synthesize liposomes. 1,2-Distearoyl-sn-glycero-3-phosphorylcholine (L-ALPHA-PHOSPHATIDYLCHOLINE, DISTEAROYL) is the lipid component in the lipid nanoparticle (LNP) system.
Targets(IC50)	Others

Solubility Information

Solubility	Ethanol: 50 mg/mL (63.28 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2656 mL	6.3279 mL	12.6558 mL
5 mM	0.2531 mL	1.2656 mL	2.5312 mL
10 mM	0.1266 mL	0.6328 mL	1.2656 mL
50 mM	0.0253 mL	0.1266 mL	0.2531 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Andrew D Miller. Delivery of RNAi therapeutics: work in progress. Expert Rev Med Devices. 2013 Nov;10(6):781-811.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481