


Cy3 diacid

Chemical Properties

CAS No. : 1034871-56-1
Formula: C₃₅H₄₅BrN₂O₄
Molecular Weight: 637.659
Appearance: no data available
Storage: keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Cy3 diacid is often used in biomolecular labeling, fluorescence imaging and other fluorescence biological analysis. They are widely used in polypeptides, proteins, nucleotides and so on.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5682 mL	7.8412 mL	15.6823 mL
5 mM	0.3136 mL	1.5682 mL	3.1365 mL
10 mM	0.1568 mL	0.7841 mL	1.5682 mL
50 mM	0.0314 mL	0.1568 mL	0.3136 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481