


7,22,25-Stigmastatrienol

Chemical Properties

CAS No. :	14485-48-4
Formula:	C29H46O
Molecular Weight:	410.686
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	7,22,25-Stigmastatrienol is a natural product
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4349 mL	12.1746 mL	24.3493 mL
5 mM	0.487 mL	2.4349 mL	4.8699 mL
10 mM	0.2435 mL	1.2175 mL	2.4349 mL
50 mM	0.0487 mL	0.2435 mL	0.487 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481