


Arenobufagin 3-hemisuberate

Chemical Properties

CAS No. :	30219-16-0
Formula:	C32H44O9
Molecular Weight:	572.69
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Arenobufagin 3-hemisuberate is a natural product
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7461 mL	8.7307 mL	17.4615 mL
5 mM	0.3492 mL	1.7461 mL	3.4923 mL
10 mM	0.1746 mL	0.8731 mL	1.7461 mL
50 mM	0.0349 mL	0.1746 mL	0.3492 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481