


Deacetylpseudolaric acid A

Chemical Properties

CAS No. : 82508-37-0
Formula: C₂₀H₂₆O₅
Molecular Weight: 346.42
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Deacetylpseudolaric acid A is a natural product from <i>Pseudolarix amabilis</i> .
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8867 mL	14.4333 mL	28.8667 mL
5 mM	0.5773 mL	2.8867 mL	5.7733 mL
10 mM	0.2887 mL	1.4433 mL	2.8867 mL
50 mM	0.0577 mL	0.2887 mL	0.5773 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481