


Delphinidin 3-rutinoside chloride

Chemical Properties

CAS No. : 15674-58-5
Formula: C₂₇H₃₁ClO₁₆
Molecular Weight: 646.98
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Delphinidin-3-O-rutinoside chloride has antioxidant activity.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5456 mL	7.7282 mL	15.4564 mL
5 mM	0.3091 mL	1.5456 mL	3.0913 mL
10 mM	0.1546 mL	0.7728 mL	1.5456 mL
50 mM	0.0309 mL	0.1546 mL	0.3091 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481