


Epivogeloside

Chemical Properties

CAS No. : 118627-52-4
Formula: C17H24O10
Molecular Weight: 388.37
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Epivogeloside is a product from <i>Lonicera japonica</i> .
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5749 mL	12.8743 mL	25.7486 mL
5 mM	0.515 mL	2.5749 mL	5.1497 mL
10 mM	0.2575 mL	1.2874 mL	2.5749 mL
50 mM	0.0515 mL	0.2575 mL	0.515 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481