


Saikosaponin G

Chemical Properties

CAS No. : 99365-19-2
Formula: C42H68O13
Molecular Weight: 780.993
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Saikosaponin G is a natural product
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2804 mL	6.4021 mL	12.8043 mL
5 mM	0.2561 mL	1.2804 mL	2.5609 mL
10 mM	0.128 mL	0.6402 mL	1.2804 mL
50 mM	0.0256 mL	0.128 mL	0.2561 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481