


2-Hydroxytetracosanoic acid ethyl ester

Chemical Properties

CAS No. : 124111-47-3
Formula: C₂₆H₅₂O₃
Molecular Weight: 412.69
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	2-Hydroxytetracosanoic acid ethyl ester is a natural product from <i>Adhatoda vasica</i>
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4231 mL	12.1156 mL	24.2313 mL
5 mM	0.4846 mL	2.4231 mL	4.8463 mL
10 mM	0.2423 mL	1.2116 mL	2.4231 mL
50 mM	0.0485 mL	0.2423 mL	0.4846 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481