


Echinoyne thiophene A

Chemical Properties

CAS No. : 64165-98-6
Formula: C₁₃H₁₀O₂S
Molecular Weight: 230.28
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Echinoyne thiophene A is a natural product from <i>Lagera pterodonta</i>
Targets(IC ₅₀)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.3425 mL	21.7127 mL	43.4254 mL
5 mM	0.8685 mL	4.3425 mL	8.6851 mL
10 mM	0.4343 mL	2.1713 mL	4.3425 mL
50 mM	0.0869 mL	0.4343 mL	0.8685 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481