


15-Isopimarene-8,18-diol

Chemical Properties

CAS No. : 73002-86-5
Formula: C₂₀H₃₄O₂
Molecular Weight: 306.48
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	15-Isopimarene-8,18-diol is a natural product for research related to life sciences. The catalog number is TN2639 and the CAS number is 73002-86-5.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2629 mL	16.3143 mL	32.6286 mL
5 mM	0.6526 mL	3.2629 mL	6.5257 mL
10 mM	0.3263 mL	1.6314 mL	3.2629 mL
50 mM	0.0653 mL	0.3263 mL	0.6526 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481