


2-Methoxyanofinic acid

Chemical Properties

CAS No. : 179457-70-6
Formula: C₁₃H₁₄O₄
Molecular Weight: 234.25
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	2-Methoxyanofinic acid is active against the plant pathogenic fungus <i>Cladosporium cucumerinum</i> .
Targets(IC50)	Antifection

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.2689 mL	21.3447 mL	42.6894 mL
5 mM	0.8538 mL	4.2689 mL	8.5379 mL
10 mM	0.4269 mL	2.1345 mL	4.2689 mL
50 mM	0.0854 mL	0.4269 mL	0.8538 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481