


20,24-Epoxy-24-methoxy-23(24-25)abeo-dammaran-3-one

Chemical Properties

CAS No. : 1020074-97-8
Formula: C₃₁H₅₂O₃
Molecular Weight: 472.754
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	20,24-Epoxy-24-methoxy-23(24-25)abeo-dammaran-3-one is used as a reference standard compound, which can be used in relevant studies in the field of life sciences.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1153 mL	10.5764 mL	21.1528 mL
5 mM	0.4231 mL	2.1153 mL	4.2306 mL
10 mM	0.2115 mL	1.0576 mL	2.1153 mL
50 mM	0.0423 mL	0.2115 mL	0.4231 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481