


3-Epikatonic acid

Chemical Properties

CAS No. : 76035-62-6
Formula: C₃₀H₄₈O₃
Molecular Weight: 456.7
Appearance: no data available
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year



Biological Description

Description	Epikatonic acid shows less active against T. (S.) cruzi trypanosome.
Targets(IC50)	Antifection

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1896 mL	10.9481 mL	21.8962 mL
5 mM	0.4379 mL	2.1896 mL	4.3792 mL
10 mM	0.219 mL	1.0948 mL	2.1896 mL
50 mM	0.0438 mL	0.219 mL	0.4379 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481