


Ageratochromene dimer

Chemical Properties

CAS No. :	17678-76-1
Formula:	C26H32O6
Molecular Weight:	440.53
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	Ageratochromene dimer is a natural product of Ageratum, Asteraceae. The catalog number is TN3358 and the CAS number is 17678-76-1. Ageratochromene dimer can be used as a reference standard.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.270 mL	11.350 mL	22.6999 mL
5 mM	0.454 mL	2.270 mL	4.540 mL
10 mM	0.227 mL	1.135 mL	2.270 mL
50 mM	0.0454 mL	0.227 mL	0.454 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481