


cis-Methylisoeugenol

Chemical Properties

CAS No. :	6380-24-1
Formula:	C ₁₁ H ₁₄ O ₂
Molecular Weight:	178.23
Appearance:	no data available
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	cis-Methylisoeugenol is a natural product of Acorus, Araceae. The catalog number is TN3666 and the CAS number is 6380-24-1. cis-Methylisoeugenol can be used as a reference standard.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.6107 mL	28.0536 mL	56.1073 mL
5 mM	1.1221 mL	5.6107 mL	11.2215 mL
10 mM	0.5611 mL	2.8054 mL	5.6107 mL
50 mM	0.1122 mL	0.5611 mL	1.1221 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 36 Washington Street, Wellesley Hills, MA 02481